

# Corsinel®

*Live as you like*

## Ostomy and Hernia Support Belt



Maximum support for hernia management and prevention

Offers discretion and relief from dragging sensation, assisting mobility

Soft, seamless and breathable for optimum comfort



Breathable



4-way stretch



Seamless




Fully Openable



Latex-free

**TYTEX®**




Corsinel Belt is designed to support and retain abdominal hernias. The belt offers maximum support to alleviate the uncomfortable feeling of a dragging hernia, making daily and strenuous activities easier to perform. Use of the belt also enhances a smooth appearance, providing discretion and encouraging participation in social life.

The belt has 4-way-stretch, which provides perfect adaption to body contour. Skin friendly material ensures no allergic reaction, redness or skin irritation in general. Seamless knitting offers optimum wearer comfort. The belt is designed for all day use.

### Indication:


The Corsinel Belt is specially designed for parastomal hernia, umbilical hernia, epigastric hernia or abdominal scar hernia. The Corsinel Belt is developed for abdominal support of the above hernia types. The support from the belt offers relief, alleviation, assists mobility and counteracts sense of gravity, thereby providing pain relief. Furthermore, the belt may help reduce possible smells and noises from a stoma. The belt may also reduce the risk of developing a hernia. The Corsinel belt is available in three heights: 15, 20 and 26 cm. For recommendation on height and use, please contact your doctor or nurse.



### Precautions:

Make sure the product is not compromising lung functions. We advise against using the belt during pregnancy, by acute abdomen, abdominal compartment syndrome, presence of ascites fluid or general discomfort/pain in the abdominal region. Please consult your doctor if you suffer from Chronic Obstructive Pulmonary Disease, if you experience unexpected discomfort, stoma or skin irritation, or if you have any questions regarding use. If urgent need for removal of product, open with scissors.

### Available colours:



### Composition:

54% Polyamide  
28% Cotton  
13% Elastane  
5% Polypropylene

### Care Handling:



Close hook-and-loop before washing. Wash with similar colours and in a washing bag. Use of chlorine based bleach or chlorine based disinfection will decrease durability and function.

### Product range:

Belt 15 cm, art. no. 3518  
Belt 20 cm, art. no. 3519  
Belt 26 cm, art. no. 3520

Measurement at the widest abdominal circumference

### Recommended size chart



SIZE	CM	IN
XS	74 - 82	29 - 32
S	78 - 90	31 - 35
M	86 - 98	34 - 39
L	94 - 106	37 - 42
XL	102 - 116	40 - 46
2XL	112 - 130	44 - 51
3XL	126 - 148	50 - 58

### 26 cm

SIZE	CM	IN
M	86 - 98	34 - 39
L	94 - 106	37 - 42
XL	102 - 116	40 - 46
2XL	112 - 130	44 - 51
3XL	126 - 148	50 - 58

NOTE: For more support go down one size and for less support go up one size.  
IMPORTANT: Always measure yourself to find the right size

Dealer: